Plasmin Inhibitor - 0020009200

Intended use

Automated chromogenic assay for the quantitative determination of Plasmin Inhibitor in human citrated plasma on IL Coagulation Systems.

Summary and principle

Plasmin Inhibitor, the major fast acting inhibitor of the fibrinolytic system, also known as Alpha₂-Antiplasmin, is an important regulator of the fibrinolytic system.¹ Congenital deficiencies are associated with hemorrhagic problems.² Decreased levels of Plasmin Inhibitor are observed in liver diseases³ and DIC.⁴ Increased levels have been reported during postoperative episodes.⁴ The Plasmin Inhibitor kit is an assay based on synthetic chromogenic substrate and on plasmin inactivation.

Plasmin Inhibitor levels in patient plasma are measured automatically on IL Coagulation Systems in two stages:

- Incubation of the plasma with the plasmin reagent in the presence of methylamine.⁵
- Quantification of the residual plasmin activity with a synthetic chromogenic substrate. The paranitroaniline released is monitored kinetically at 405 nm and is inversely proportional to the Plasmin Inhibitor level in the test sample.

Composition

The **Plasmin Inhibitor** kit consists of:

[B] **Buffer** (Cat. No. 0020009230): 2 x 9 mL vials of a concentrated buffer solution containing sodium chloride, methylamine and surfactant.

[E] **Plasmin reagent** (Cat. No. 0020009220): 2 x 2.5 mL vials of a lyophilized preparation containing human plasmin (2.5 nkat/vial), buffer, human serum albumin, stabilizers and bulking agent.

[S] **Chromogenic substrate** (Cat. No. 0020009210): 1 x 4 mL vial of the lyophilized chromogenic substrate S-2403, pyroGlu-Phe-Lys-pNA.HCl (8 mg/vial) and bulking agent.

PRECAUTIONS AND WARNINGS:

The material in this product was tested with FDA cleared methods and found nonreactive for Hepatitis B surface Antigen (HBsAg), Anti-HCV and HIV antibodies. Handle as if potentially infectious.⁶

Hazard class: **None**
Risk phrases: **None**
Safety phrases: **None**

This product is For *in vitro* Diagnostic Use.

Preparation

Buffer: Dilute the necessary quantity of the concentrated buffer 1:10 (1+9) with CLSI Type CLR water or equivalent.⁷ Mix before use.

Plasmin reagent: Dissolve the contents of each vial with 2.5 mL of Diluted Buffer. Replace the stopper and swirl gently. Make sure of the complete reconstitution of the product. Keep the reagent at 15-25°C for 30 minutes and mix gently before use.

Chromogenic substrate: Dissolve the vial contents with 4 mL of CLSI Type CLR water or equivalent.⁷ Replace the stopper and swirl gently. Make sure of the complete reconstitution of product. Keep the substrate at 15-25°C for 30 minutes and mix gently before use.

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Reagent storage and stability

Unopened reagents are stable until the expiration date printed on the vial when stored at 2-8°C.

Buffer - Opened reagent should be kept at 2-8°C in the original vial.

Diluted Buffer - Stability after dilution: 24 hours at 15°C.

Plasmin reagent - Stability after reconstitution: 5 days at 15°C and 2-8°C, 3 months at -20°C in the original vial, 5 days at 15°C on the ACL ELITE[®]/ ELITE PRO/8/9/10000 or 24 hours at 15°C on the ACL Futura/ACL Advance and ACL TOP[®] Family.

Chromogenic substrate - Stability after reconstitution: 5 days at 15°C and 2-8°C, 3 months at -20°C in the original vial, 5 days at 15-25°C on the ACL ELITE[®]/ ELITE PRO/8/9/10000 or 24 hours at 15°C on the ACL Futura/ACL Advance and ACL TOP Family.

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial.

Instrument/test procedures

Refer to the appropriate IL instrument's Operator's Manual and/or Application Manual for the complete assay procedure instructions.

Specimen collection and preparation

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Refer to CLSI Document H21-A5 for further instructions on specimen collection, handling and storage.⁸

Additional reagents and control plasmas

The following are not supplied with the kit and must be purchased separately.

	Americas and Pacific Rim	Europe
	Cat. No.	Cat. No.
Calibration plasma	0020003700	0020003700
Normal Control	0020003120	0020003110
Special Test Control Level 1	0020011000	0020011000
Special Test Control Level 2	0020012000	0020012000

Quality control

Normal and abnormal controls are recommended for a complete quality control program.⁹ Normal Control and Special Test Control Level 1 and 2 are designed for this program. Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance with good laboratory practice. Refer to the instrument's Operator's Manual for additional information. Refer to Westgard *et al* for identification and resolution of out-of-control situations.¹⁰

Results

Plasmin Inhibitor results are reported in % activity. Refer to the instrument's Operator's Manual for additional information.

Limitations/interfering substances

Plasmin Inhibitor results on the ACL, ACL ELITE/ELITE PRO/ 8/9/10000 and ACL Futura/ ACL Advance Systems are not affected by heparin (UF heparin or LMW heparin) up to 2 U/mL, α₂-macroglobulin up to 7 mg/mL, hemoglobin up to 200 mg/dL, bilirubin up to 20 mg/dL and triglycerides up to 1000 mg/dL.

Plasmin Inhibitor results on the ACL TOP Family are not affected by heparin (UF heparin or LMW heparin) up to 2 U/mL, α₂-macroglobulin up to 3.5 mg/mL within the Plasmin Inhibitor normal range, hemoglobin up to 500 mg/dL, bilirubin up to 24 mg/dL and triglycerides up to 2320 mg/dL.

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ENGLISH - Insert revision 05/2013

Expected values

A normal range study was performed using Plasmin Inhibitor kit.

System	N	Range (units)
ACL Family	50	89,3 - 112,1 (% activity)
ACL Futura/ACL Advance	60	96,8 - 118,8 (% activity)
ACL TOP Family	124	98,0 - 122,0 (% activity)

Ranges were calculated as recommended by the International Federation of Clinical Chemistry (IFCC).¹¹ These results were obtained using a specific lot of reagent. Due to many variables which may affect results, each laboratory should verify its own normal range.

Performance characteristics

Precision:

Within run and total (run to run and day to day) precision was assessed over multiple runs using both normal and abnormal samples.

ACL Family	Mean (% activity)	CV % (Within run)	CV % (Total)
Normal Control	102	2.0	2.5
Abnormal Control	46.6	2.7	3.9

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- Collen D, Wiman B. Fast-acting Plasmin Inhibitor in Human Plasma, Blood. 1978; 51: 536-569.
- Booth NA. The Natural Inhibitors of Fibrinolysis, Haemostasis and Thrombosis. Editors, Bloom AL, Forbes CH, Thomas DP, Tuddenham EGD. Edition 3. 1994; 1: 699-717.
- Aoki N, Yamanaka T. The α₂-Plasmin Inhibitor Levels in Liver Diseases, Clin. Chem. Acta. 1978; 84: 99-105.
- Matsuda M, Wakabayashi K, Aoki N, Morioka Y. α₂-Plasmin Inhibitor is among Acute Phase Reactants, Thromb. Res. 1980; 17: 527-532.
- Billing Clason S, Meijer P, Kluff C, Ersdal E. A Novel Method COAMATIC® Plasmin Inhibitor, for the Specific Determination of Plasmin Inhibitor Activity in Plasma, Fibrinolysis. 1996; 10. Suppl. 2: 165-167.
- Richmond JY, McKinney RW eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, 4th Edition, 1999.

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